



Reach for the stars with StarzTipz®

Colourful – Versatile – Customisable – Reduces Cross-contamination



We will help you identify your syringe* and supply a conversion kit (if required) to make the switch to StazTipz as easy as possible.

*see QR code on the back page

			
Conversion Kit Cefla Anthos F3/F6 QRNT1030	Conversion Kit Cefla SP QRNT1035	Conversion Kit Faro SM03 QRNT1100	Conversion Kit Kavo Estetica QRNT1092
			
Conversion Kit Kavo K4 QRNT1090	Conversion Kit Luzzani Minibright QRNT1055	Conversion Kit Luzzani Minilight QRNT1050	Conversion Kit Luzzani Minimate QRNT1051
			The item numbers only refer to StarzTipz® adapters
Conversion Kit Quick Release ADEC/UNIC QRNT1000	Conversion Kit Sirona Sprayvit 4000 QRNT1071	Conversion Kit Sirona Sprayvit C8 QRNT1070	

Find answers to the most frequently asked questions about StarzTipz®. Our aim is to provide you with all important information quickly and clearly. Should you have any further questions, we will be happy to answer them personally.

01. What is the advantage of our “6 Air Channel“ design?

By having 6 small air channels, the dental air is more concentrated and dispensed through each of the 6 channels, giving you a more focused spray through each channel. That equates to a more powerful dental air spray.

02. What is special about the HP tip?


With its length of 65 mm and double chambers, the High-Pressure (HP) tip delivers a more powerful and focussed spray and is ideal for cleaning deep periodontal pockets. It is 30 mm shorter than the original tip and is ideal for use with young patients and for working close to the treatment site.

03. Does this tip work on my syringe?

We have a range of converters to fit many syringes.

04. Can I really be sure that StarzTipz® guarantees an exact separation of air and water?

If your equipment is in good working order, we guarantee 100% separation of air and water. The StarzTipz® attachment sits firmly on the syringe thanks to the O-ring lock. In addition, StarzTipz® translucent outer core allows you to see if water has entered into the air chamber. 99% of leakage problems are due to worn O-rings, incorrect fitting of the tip or from trying to fit the tip with an incompatible adapter.



StarzTipz® Key Features

- ★ Immediate separation of air and water
- ★ Flexible plastic, retains its shape
- ★ Worldwide compatibility through matching adapter kits
- ★ Slim design for optimum visibility into the oral cavity
- ★ Reduced risk of cross-contamination

Switching to StarzTipz® is easier than you think!

Select the model of your air/water syringe and find the right adapter for you and the corresponding O-ring kit to convert your syringe system to StarzTipz®.



Switch to
StarzTipz® now!